

Electrical Bed Series

This Electric ICU Bed, Model B6001 series has been developed throughout the knowledge, skilled personnel and advance technology in this field. We are dedicated to develop the most "Advance Ultimate innovation" products by a nonstop research to ensure the excellent quality of engineering work.

B6001 series is a ICU/CCU electrical hospital bed with a latest technologies which provides a full force in nursing care.



Description

The beds to facility significantly eases the physical strain on nursing staff helping them perform routine nursing activities with the minimum of effort. Routine nursing activities often mean the physical handling of the patient several times a day, often in non-ergonomic positions.

PRODUCT SPECIFICATION & PERFORMANCE




1. Ventilated 4-section mattress support base.
2. Removable and detachable PP-moulded head and foot board foot board.
3. Central locking or castor with individual lock.
4. Provision for IV Pole attachment x 4/6 nos

5. Rotation bumpers or built in bumper bu at all 4 corners
6. Complete with or without X-ray cassette holder at the backrest section
7. Complete with mattress or pressure reducing mattress with hospital grade mattress cover (anti bacterial, anti-static, water proof and fire retardant).



1. Electric motorized bed with 4 motors 5 function
2. With or without auto regression for the backrest and legrest ; if with auto regression is to avoid shear and pressure point on patient.
3. 2 pcs of PE side rail with embedded control button (customarization with upon customer request)
4. HI-LO, backrest ($75^{\circ} \pm 5^{\circ}$) & footrest ($35^{\circ} \pm 5^{\circ}$) elevation operated by Electric Linear Actuator. Trendeleburg ($12^{\circ} \pm 2^{\circ}$) and reverse trendelenburg ($12^{\circ} \pm 2^{\circ}$) by electric motorized
5. With or without back-up battery, working up to 20-25 circle.
6. With or without care giver function (nurse control and nurse lock.)
7. With or without cardiac chair position.



Excellent in Medical Technology

Advantages of B6001-8VV / 8VVW



Double Auto-Regression 10+9 cm-Latest Technology for model with this function only.

When repositioning the backrest and thigh rest, makes the back plate retreat 10cm and the thigh plate retreats 6cm so to minimises the pressure and shear forces in the patient's pelvic region.

This results in a significant reduction in the development of pressure ulcers.



Degree Indicator Embedded to side rail



Bed Extension 200mm



Pressure Reducing Mattress



Central Locking Wheel with direction lock



Stainless steel IV Pole-4 hooks



Available for weighing scale model only
Maximum weighing range 500kg/100lb,
alarm accuracy $\pm 0.1\text{kg} / \pm 0.2\text{lb}$



Manual CPR emergency reset handle
to immediately flatten the bed back
part; located on both sides of the bed.



Safety
Follow IEC60601 safety standard one
button emergency stop switch can
immediately stop all going movement,
for applicable model.

**Excellent in Medical
Technology**

Electrical Bed Series

Model :
B6001

Description :
ICU/HDU/Ward 5
Function Bed

Sample :



Model :
B6001-8VV

Description :
ICU/HDU/Ward 5
Function Bed

Sample :



Model :
B6001-8VVW

Description :
ICU/HDU/Ward 5
Function Bed

Remarks :
With weighing scale



Sample :



Model :
B6001-8VvV

Description :
ICU/HDU/Ward 5
Function Bed

Remarks :
With vibration
function



Sample :



Electrical Bed Series

Model :
B6001-8YT

Description :
ICU/HDU/Ward 5
Function Bed

Sample :



Model :
B6001-8YTW

Description :
ICU/HDU/Ward 5
Function Bed

Remarks :
With weighing scale



Sample :



Model :
B6001-8YY

Description :
ICU/HDU/Ward 5
Function Bed

Remarks :
Removable plastic covers blow
moulded to minimize the slots
where unwanted infections
may hide



Sample :



Model :
B6001-FLSR

Description :
ICU/HDU/Ward 5
Function Bed
– Full Length epoxyside rail

Remarks :
Metal siderail
– full length

Sample :



Electrical Bed Series

Model :
B6001-6Y

Description :
ICU/HDU/Ward 5
Function Bed

Sample :



Model :
B6001-8f

Description :
ICU/HDU/Ward 5
Function Bed

Remarks :

Nurse lock embedded at
footboard with LCD display
screen

Sample :



Model :
B6001-8YYW

Description :
ICU/HDU/Ward 5
Function Bed

Remarks :

With weighing scale

Sample :



Model :
B6001-SD

Description :
ICU/HDU/Ward 5
Function Bed

Remarks :

Control button
embedded to siderail

Sample :



Electrical Bed Series

Model :
B6001-WS

Remarks :
With built in weighing scale

Sample :

Description :
ICU/HDU/Ward 5
Function Bed



Model :
B6003

Remarks :
Bariatric use

Sample :

Description :
Electric Bed 5
function



Model :
B6007

Remarks :
Bariatric use

Sample :

Description :
Electric Bed 5
function



Model :
B5001-OF

Remarks :

Sample :

Description :
Electric Bed/Ward
3 Function



Electrical Bed Series

Model :
B5001-BS

Remarks :

Sample :

Description :
Electric Bed/ Ward
3 motors 4 function



Model :
B5001-SC

Remarks :

Sample :

Description :
Electric Bed/ Ward 3
motors 4 function



Model :
B5001-6VK

Remarks :

Sample :

Description :
Electric Bed/ Ward 3
motors 4 function



Model :
B5001-6VV

Remarks :
PP siderail with
extra height

Sample :

Description :



Electrical Bed Series

Model :
B5001-TCPR

Remarks :
Full length aluminum
siderail with the filler
for ¼ gap

Sample :

Description :



Model :
B5001-CPR

Remarks :

Sample :

Description :



Model :
B5001-FLSR

Remarks :
12 strong J Bar siderail
cover up 90% of bed

Sample :

Description :
Electric Bed/
Ward 3 Function



Model :
B5001-PE

Remarks :
PE siderail

Sample :

Description :
Electric Bed/ Ward 3
Function



Electrical Bed Series

Model :
B5001-VS

Remarks :

Sample :

Description :



Model :
B5002

Remarks :

Sample :

Description :
Electric Bed –
Ultra Low



Model :
B5002-NL

Remarks :

Sample :

Description :
Electric Bed –
Ultra Low



Model :
B5000

Remarks :

Sample :

Description :
Electric Bed/ Ward 4
function



Electrical Bed Series

Model :
B5000-6YT

Remarks :

Sample :

Description :



Optional Accessories

Model : RX100	Remarks : Shoe Compartment	Sample : 
Description : Hospital Bedside Rack		

Model : RX100-1	Remarks :	Sample : 
Description : Hospital Bedside Rack		

Model : RX100-1SL	Remarks : Smart Locker	Sample : 
Description : Hospital Bedside Rack	 	

Model : RX101	Remarks : Mild Steel with epoxy coated	Sample : 
Description : Hospital Bedside Rack		

Optional Accessories

Model :
RX101-G

Remarks :
3 ways guard rail

Sample :

Description :
Hospital Bedside
Rack



Model :
RX101-TGR

Remarks :
3 sided top safety
guardrail

Sample :

Description :
Hospital Bedside
Rack



Model :
RX103

Remarks :

Sample :

Description :
Hospital Bedside
Rack



Model :
RX104

Remarks :
Laminated chipboard
constructed

Sample :

Description :
Hospital Bedside
Rack



Optional Accessories

Model :
RX104-SL

Remarks :
Smart Locker

Sample :

Description :
Hospital Bedside
Rack



Model :
RX201

Remarks :
Adjusted by classic
coil spring for Hi-Lo
adjustment

Sample :

Description :
Hospital Cardiac
Table



Model :
RX202

Remarks :
Operated by 2 finger
for Hi-Lo adjustment

Sample :

Description :
Hospital Cardiac
Table



Model :
RX202-C

Remarks :
Come with drawer
and mirror underneath

Sample :

Description :
Hospital Cardiac
Table



Optional Accessories

Model :
RX202-D

Remarks :
Come with drawer
and mirror underneath

Sample :

Description :
Hospital Cardiac
Table



Model :
RX202-T

Remarks :

Sample :

Description :
Hospital Cardiac
Table



Model :
RX203

Remarks :
Gold plated sheet

Sample :

Description :
Hospital Cardiac
Table



Model :
RX205

Remarks :
C shape leg which is
allow the table to go
in to dialysis chair

Sample :

Description :
Hospital Cardiac
Table



Optional Accessories

Model :

RXHP1

Remarks :

Sample :

Description :

Pillow Adult



Model :

RX202-T

Remarks :

Sample :

Description :

Pillow Adult



Model :

P-Slide

Remarks :

Sample :

Description :

Patslide



Model :

CHPU01-1

Remarks :

Sample :

Description :

Patient Chair



Optional Accessories

Model :
CHPU02

Remarks :

Sample :

Description :
Patient Chair




Model :
CHPU03-B

Remarks :
Bariatric use

Sample :

Description :
Patient Chair




Model :
CHPU03

Remarks :

Sample :

Description :
Patient Chair



Model :
CHR04-2

Remarks :

Sample :

Description :
Recylining /
Hemodialysis chair



Optional Accessories

Model :
CHR06

Remarks :

Sample :

Description :
Recylining /
Hemodialysis chair



Model :
CHR07

Remarks :
Removeable ABS meal
table and able to
secure at side

Sample :

Description :
Geriatric Chair /
Recylining Chair



Model :
CHR09

Remarks :

Sample :

Description :
Electric Recylining
Chair



Model :
CHR11

Remarks :

Sample :

Description :
Reclining Chair



Optional Accessories

Model :
CHSP02

Remarks :

Sample :

Description :
Accompany Chair /
Sleeper Chair



Model :
CHSP02-B

Remarks :
Removeable ABS meal
table and able to
secure at side

Sample :

Description :
Accompany Chair /
Sleeper Chair



Model :
CHSP05

Remarks :

Sample :

Description :
Accompany Chair /
Sleeper Chair



Model :
CHSP05-EW

Remarks :

Sample :

Description :
Accompany Chair /
Sleeper Chair



Optional Accessories

Model :
CHWC01

Remarks :
Complete with Iv
pole hook

Sample :



Description :
Wheelchair

Model :
CHWC02

Remarks :
Complete with Iv pole
hook

Sample :



Description :
Wheelchair

