


Linen Cart

Model :
CLBS03

Remarks :

Sample :

Description :
Hospital
Sponging Trolley



Model :
CLCD32

Remarks :

Sample :

Description :
Clean & Dirty
Linen Cart



Model :
CLD01

Remarks :

Sample :

Available:
CLD02
CLD03

Description :
Linen Distribution



Model :
CLDB2

Remarks :

Sample :

Description :
Double Linen Bag



Linen Cart

Model :
CLS01

Remarks :

Sample :

Description :
Hospital
Linen Cart



Model :
CLSB1

Remarks :

Sample :

Description :
Single Linen Bag



Model :
SLT01

Remarks :

Sample :

Description :
Soiled
Linen Trolley



Model :
SLT02

Remarks :

Sample :

Description :
Double Soiled
Linen Trolley



Linen Cart

Model :
TSP02

Remarks :

Sample :

Description :
Sponging Trolley



CSSD

Model :
C9050

Remarks :

Sample :

Description :
Distribution Cart



Model :
PSTD2

Remarks :

Sample :

Description :
Paper Stand



Model :
CSSDWT1260

Remarks :
Available :
CSSDWT1575
CSSDWT25100

Sample :

Description :
CSSD
Working Table



Model :
PT15075

Remarks :

Sample :

Description :
Preparation Table



CSSD

Model :
PTH500

Description :
Pass
Through Hatch

Remarks :
Available:
PTH600
PTH700
PTH800

Sample :



Model :
WT6030-D

Description :
Working Table
with Drawer

Remarks :

Sample :



Model :
UTRSH1

Description :
Sloop Hopper
with sink
Table top with
marrine

Remarks :
Complete with
working table

Sample :



ICU Chart Table (COW) DDA

Model :
DDA180

Remarks :
Dimension (mm) :
715W x 180D x 760H

Sample :

Description :
Hospital
DDA Cabinet

Battery operated
buzzer



Model :
CTB100

Remarks :
Mild Steel sheet for
the framework

Sample :

Description :
ICU Chart Table

Front table top
with higher lips



Model :
LTS02

Remarks :
Dimension (mm) :
520W x 420D x 860-1050H

Sample :

Description :
Computer on
Wheel (COW)

Come with sharp
bin holder



Model :
PBS01

Remarks :
Complete with 3
layer of drawer

Sample :

Description :
Phebletomy Cart



Mortuary Items

Model :
MOAT-A02

Description :
Autopsy Table

Remarks :

Specifaction :
Material: Stainless steel 304
Product size (mm): L2600 x W800 x H800
Product size (H): L102.3 x W31.5 x H31.5
Table Height Adjustable (mm): 720-920 by motor
The mobile spray gun length(mm): 2000
Power: 560W
Horsepower: 1.0hp
Speed: 4200r/Min
Capacity: 1200ML
Emergency stop button: one (control main power off)
Handheld remote: yes
Voltage: 220-240V/50Hz; Electricity 1.5A
Load capacity: 200KG
Package size(mm): 2630 x 830 x 830

Sample :



Model :
BPTSS01

Description :
Body
Preparation Table

Remarks :
For body preparation,
Material with 304
grade stainless steel
and load capacity to 200kg

Dimension (mm):
2010L x 710W x 750

Sample :




Model :
MOEL1530

Description :
Lift Trolley

Remarks :

Material	Stainless steel 304
Operation	Fully electric (Battery powered)
Operation method	Remote key
Scissors	2
Adjustable lift height(mm):	450 to 2000
Dimensions lowest position W/D/H:	2060*670*450
Dimensions highest position W/D/H:	2060*670*2000
Number of trays	1 Piece
Number of handles on tray	One on each side. (Customization available)
Body tray size(mm)	1940*580
Thickness of stainless steel	1.0mm
Carrying capacity	200kg
Wheel	4pcs (2pcs locked)
Single person operation	Yes
Battery	12v
Charging time	4-8 hours
Package	Plywood
Packing Size(mm); Weight	2230*710*560; 138kg
Certification	CE
Customization	Yes

Sample :



Model :
MOAT2-DT004

Description :
Mortuary Trolley

Remarks :
Dimension (mm):
2500L x 800W x 850H

Include:
Sink 1 pcs
Socket safety
Water tap
Water spray gun

Sample :



Mortuary Items

Model :
MOAT3-DT003

Remarks :
Dimension (mm):
2400L x 700W x 850H

Description :
Mortuary Trolley

Include:
Drain pipe 1 pcs
Castor 4 pcs

Sample :



Racking System


Model :
RKBL1218-4

Description :
Rack Boltless

Remarks :
Multipurpose rack with
4 pcs fibreboard shelve

Dimension (mm) :
1200W x 1800H x 600D

Sample :



Model :
SSR05

Description :
Multi Purpose
Rack

Remarks :
5 shelves

Dimension (mm):
1410L x 520W x 1600H

Sample :



Model :
RKSWM4

Description :
Multi Rack
Wire Mesh

Remarks :
Dimension (mm) :
1570L x 610D x 1813H

Sample :



Utilities Rooms

Model :
CWM01

Remarks :
CSSD rack with
wire mesh shelving

Sample :



Description :
CSSD Rack

Model :
OTBD1

Remarks :
Brus dispenser – wall
mounted

Sample :




Description :
Brush Dispenser

Dimension (mm) :
700W x 125D x 610H

Model :
UTRB1

Remarks :
Dimension (mm) :
290L x 96W x 90H

Sample :



Description :
Catheters Basket

Model :
UTRSCS1

Remarks :
Fully stainless steel
constructed,
304 grade

Sample :



Description :
Cleaner sink

Utilities Rooms

Model :
UTRSH1

Description :
Sloop Hopper

Remarks :
Wall mounted /
floor standing

come with 1 slopper
and 1 sink
underhearth with
storage shelve

Sample :



Model :
UTRW01

Description :
Mop Holder
with Shelve

Remarks :
Wall mounted

come with 3 rubber
holder for mop with
spring lock and
4 utility hook



Sample :



X-Ray Department

Model :
XRVL101-SB

Description :
X-Ray Viewer
LED with single
viewer

Remarks :
Available for
XRVL102-SB
XRVL103-SB
XRVL104-SB

Sample :



Model :
MXVB1

Description :
X-Ray Viewer
with single viewer

Remarks :
Available for
MXVB2
MXVB3
MXVB4

Sample :




Model :
XRAH01

Description :
Apron Hanger

Remarks :
Available for wall
mounted with
model XRAHW1

Sample :



Model :
XRYLS1

Description :
Lead Screen

Remarks :
Available for
XRYLS2



Sample :

