

Avunncare Bed Series

Model :
AVÜNNCARE-eB

Description :
AVÜNNCARE-eB

is a four-section electric profiling patient bed with lateral tilt is ideal for the care of immobile patients in an array of hospital setting.

Back section adjustment:

- Anti Clamping Function
- Backrest 30° stop function
- Backrest buffer function
- Auto regress function
- X-ray translucent platform
- Auto contour function
- Mattress function
- Weighing function
- Alarm function
- Night light function
- One button reset function
- Emergency stop function

Sample :



Function:

- Backrest max upward angle: 75°
- Footrest max upward angle: 40°
- Height is adjustable: 480 – 880mm
- Max Reverse – Trendelenburg: 0 – 12°
- Max Trendelenburg : 0 – 12°
- Lateral Tilt Left : 0 – 15°
- Lateral Tilt Right : 0 – 15°
- Electric Bed Extension: 150mm
- Auto Counter Function
- Auto Regress Function
- One Button Cardiac Chair Position Function
- One Button Electric CPR Function
- Manual CPR Release Accessible On Both Sides
- One Bottom Chair Position Function



Lateral Tilt left & right



Leg section up & down



Avunncare Bed Series

Model :
AVÜNNCARE

(Lateral Tilt ICU Bed with manual castor movement)

Description :
AVÜNNCARE

is an advance ICU/CCU electrical hospital bed with a lateral tilt which provides a full force in nursing care

7 advantage of AVÜNNCARE :

- Reduce the physical strain on nurses when handling a patient
- Eliminates the risk of occupational strain and physical disability
- Save time required for routine nursing activities
- Makes prevention of pressure sores more effective
- Helps prevent complication in patient with lymphoedmas and facilities post operative drains
- Helps prevent some negative signs of the immobilization syndrome
- Helps increase the quality of the life of long – term or permanently immobile patient.

Sample :

Remarks :

External dimensions 2180 × 1040mm
 Mattress dimensions 1980*900*100mm
 Mattress platform 1980 x 900mm
 height adjustment 500–800mm
 Lateral tilt 20°
 Trendelenburg / Anti Trendelenburg ± 15°
 Backrest tilt 0°–70°
 Thigh rest tilt 0°– 30°
 Safety work load 230 kg
 Weighing Scale Capacity: 300kg
 Weighing Scale Tolerant 100mg
 Degree indicator Tolerance :±2°



Avunncare Bed Series

Model : AVÜNNCARE

(Lateral Tilt ICU Bed with electric castor movement)

Description :

Electric Bed with electric Motorized wheel is more advance model with electric castor that uses electric transportation castor to minimize number of staff nurse while transporting patient and lesser workload. It is with 2 hands operation on the ergonomic handles with right hand handle movement control using fingertip control button.



The electrical castor mounted at middle of the chassis that acts as 5th wheel of the bed to make the bed with single operator during the transportation, easy maneuver and hassle free IP66 electric castor.

The higher torque will allow the castor to handle all kinds of rough surfaces.

Remarks :

External dimensions 2180 × 1040mm
Mattress dimensions 1980*900*100mm
Mattress platform 1980 x 900mm
height adjustment 500-800mm
Lateral tilt 20°
Trendelenburg / Anti Trendelenburg ± 15°
Backrest tilt 0°-70°
Thigh rest tilt 0°- 30°
Safety work load 230 kg
Weighing Scale Capacity: 300kg
Weighing Scale Tolerant 100mg
Degree indicator Tolerance :±2°

Sample :

