

Hospital Manual Crank Bed Series

Model :
B2001, Single Fowler
(Backrest adjustable
by manual cranking)

Description :
Hospital Bed
– Manual Crank Fixed
Height Bed

Remarks :
Specification:
Material: Mild Steel
Framework
Treatment: Zinc
Phosphate Washing
Finishing: Epoxy Powder
Coating
Dimension (mm) :
2160L x 900W x 500H
Nett Weight: 60kg
Load Capacity: 250kg

Sample :



Model :
B2002, Double Fowler
(Backrest adjustable
by manual Cranking)

Description :
Hospital Bed –
Manual Crank
Fixed Height Bed

Remarks :
Specification:
Material: Mild Steel
Framework
Treatment: Zinc
Phosphate Washing
Finishing: Epoxy Powder
Coating
Dimension (mm) :
2160L x 900W x 500H
Nett Weight: 60kg
Load Capacity: 250kg

Sample :



Model :
B3001, Single Fowler
(Hi-Lo, Backrest only
by manual cranking)

Description :
Hospital Bed –
Manual Crank Fixed

Remarks :
Specification:
Material: Mild Steel
Framework
Treatment :
Zinc Phosphate
Washing
Finishing: Epoxy Powde
Coating
Dimension (mm) :
2160L x 950W x 420
-740 ± 20mm
Nett Weight: 80kg
Load Capacity: 200kg

Sample :



Hospital Manual Crank Bed Series

Model :
B3002, Double Fowler
(Hi-Lo, Backrest and
legrest by manual
Cranking)

Description :
Hospital Bed –
Manual Crank Fixed

Remarks :
Specification:
Material: Mild Steel
Framework
Treatment: Zinc
Phosphate Washing
Finishing: Epoxy Powder
Coating
Dimension (mm) :
2160L x 950W x 420
-740 ± 20mm
Nett Weight: 80kg
Load Capacity: 200kg

Sample :



Model :
B3002 – 3KK

Description :
Hospital Manual
Crank Bed

Remarks :
Specification:
External size (mm)
:2160L x 980W x 400
-710 ±10mm
Mattress Platform (mm):
1940L x 830 ±10W
Backrest Adjustment :
0-85°±10°
Knee rest Adjustment :
0-45°±10°
Height Adjustment : 400-
710mm±10mm

Sample :



Model :
B3002-3VV

Description :
Hospital Manual
Crank Bed

Remarks :
Specification:
External size (mm) :
2130L x 1025W
Mattress Platform
(mm): 1925L x 900W
Safe Working Load :
250kg
Height adjustment :
430-740mm
Back rest adjustment :
0 - 73° (±5°)
Auto regression: 100mm
PE Siderail

Sample :



Hospital Manual Crank Bed Series

Model :
B3002-3AK

Description :
Hospital Manual
Crank Bed

Remarks :
Specification:
External size (mm) :
2110L x 960W x 400
-710H
Safe working Load :
200kg
Height adjustment :
400-710
Back rest adjustment :
0 - 75° (±10°)

Sample :



Model :
B4001

Description :
Hospital Bed –
Hydraulic

Remarks :
Specification :
Material : Mild Steel
Framework
Treatment : Zinc –
Phosphate Washing
Finishing : Epoxy
Powder Coating
Dimension (mm) :
2130L x 975W x 475-780H
Nett Weight : 160kg
Load Capacity : 250kg

Sample :



Model :
B4002

Description :
Hospital Bed –
Hydraulic

Remarks :
Specification :
Material : Mild Steel
Framework
Treatment : Zinc
Phosphate Washing
Finishing : Epoxy Powder
Coating
Dimension (mm) : 2150L x
1000W x 450-700
Nett weight: 160kg
Load Capacity: 250kg

Sample :



Hospital Manual Crank Bed Series

Model :
Siderail

Remarks :

Sample :

Description :
90% Collapsible
siderail



Model :
Siderail

Remarks :

Sample :

Description :
Epoxy siderail –
full length sliding
up and down

