

Women & Children

Model :
LDR11

Description :
Hospital Labour Bed
– Fixed Height

Remarks :
Specification :

Material: Mild Steel Framework

Treatment : Zinc Phosphate Washing

Finishing : Epoxy Powder Coating

Dimension (mm): 1830L x 760W x 750H

Load Capacity: 180kg

Accessories :

- 1 unit stainless steel, double hooked adjustable IV pole
- 1 unit stainless steel soiled basin (360mm) mounted on a retractable basket.
- 1 pair adjustable padded lithotomy stirrups

Sample :



Women & Children

Model :
LDR20

Description :
Hospital Labour Bed
– Hydraulic

Remarks :
Specification :

Material :

Mild Steel framework

Treatment : Zinc Phosphate Washing

Finishing : Epoxy Powder Coating

Dimension (mm) : 1340 – 1920L x 1030W
x 560–830H

Load Capacity: 250kg

Accessories :

- Stianless Steel adjustable hand grip
- Stainless steel adjustable legrest
- Knee Crutch
- Squatting support

Sample :



Women & Children

Model :
LDR31

Description :
Hospital Labour

Remarks :
Specification :

Material :

Dimension (mm) : 2025L x 1045W
x 525-960 ± 20mm

Backrest : 70°±10°

Legrest: Retractable underneath

Tredelenburg: 21°±2°

Reverse Trendeleburg: 16°±2°

Corner Bumper: yes, 4 corner

Accessories :

- Screen holder during labor pain
- knee crutch
- Footing support
- IV pole with 4 hooks

Sample :



Women & Children

Model :
LDR32

Sample :

Description :
Hospital Labour Bed

Remarks :

Specification :

Material :

Dimension (mm) :

1950 ± 20L x 860 ± 15W

Back section turning : >63°

Seat section turning: >15°

Foot rest rotate out: >90°

Foot rest turn up: >90°

Siderail lifting (mm): 400±20

Load capacity: 280kg

Accessories :

- 2 pcs legrest
- 4 pcs clamp
- 2 pcs Grab handle
- 2pcs Foot rest
- 1pcs Filth basin , Cushion, IV stand



Women & Children

Model :
NEBBT1

Sample :

Description :
Hospital Baby Bath Tub



Remarks :

- Framework : Stainless steel constructed
- Come with 1 cabinet and 1 open shelf
- Table top come with 1 plastic bath tub and side with baby changing cushion PVC padded

Model :
NEBC10

Sample :

Description :
Hospital Baby Bassinet



Remarks :

- Framework : Stainless steel with hard epoxy powder coated
- come with 1 open shelf
- Infant cot : made of clear acrylic
- come with PVC padded cushion

Women & Children

Model :
NEBC11

Sample :

Description :
Hospital Baby Bassinet

Remarks :

- Framework : Stainless steel constructed with hard epoxy powder coated
- come with 1 cabinet and 1 open shelf
- come with PVC padded cushion



Model :
NEBC30B

Sample :

Description :
Hospital Baby Bassinet

Remarks :

- Height adjustable by gas spring, easy for operate
- Transparent ABS baby basin
- Load Capacity: 25kg
- Basin lifting: 565-690mm
- Basin left /right tilting : 0~10°

NEBC30B – Blue

NEBC30G – Green

NEBC30P – Pink



Women & Children

Model :
NEBC23

Sample :

Description :
Hospital Baby Bassinet

Remarks :

- With safety lock and able to tilt up to $+10^{\circ}$ / -10° ☒- Come with 1 drawer and 1 cabinet
 - come with 2 pushing handle
- PVC padded cushion
Dimension (mm) : 905L x 500W x 950H

